

CATALOG



The JOINSTAR COVID-19 Antigen Layman's Rapid Test is approved for self-administration by laymen and is used for the qualitative detection of SARS-CoV-2 antigens. A simple and painless sample collection from the anterior nasal region

makes the handling very easy and the result can be read after 15 minutes.

Short Description:

- Joinstar COVID-19 ANTIGEN RAPID TEST
- Easy handling, per rod in the front nose area
- Packs of 1 and 5
- Layman test for self application
- Test officially approved and suitable for laypersons by the BfArM.

Intended use:

The JOINSTAR COVID-19 Antigen Layman's Rapid Test is approved for self-administration by laymen and is used for the qualitative detection of SARS-CoV-2 antigens. A simple and painless sample collection from the anterior nasal region makes it foolproof to use and the result can be read in 15 minutes.

Components:

- SARS-CoV-2 antigen test cassette
- Extraction solution in tube
- Instructions for use
- Collection Swab

Properties:

- Sensitivity = 96.10%
- Sensitivity = 96.10%
- Easy to use - Fast and reliable test results in only 15 minutes
- BfArM listed with special approval
- Manufacturer: Joinstar Biomedical Technology Co. LTD, Hangzhou, Zhejiang, China.


Further notes

- The kit should be stored at room temperature 2 to 30°C maximum.
- Sensitivity of the test decreases after the first 5 days after the onset of symptoms
- Test can only tell if antigen is present, not in what concentration
- Test result depends on quality of specimen collection, processing, transport and storage. Any errors may cause inaccurate test results.
- People under 18 years of age should be assisted by an adult
- Intended for in vitro diagnostic use only
- Please be sure to follow the manufacturer's information within the

instructions for use and on the product packaging

[Read More](#)

Category: [Corona Self Tests](#)



The Hotgen Coronavirus Antigen Test is imported by us directly from the manufacturer.

It is for single use and is used for the determination of SARS-CoV-2 antigen using a human nasal swab.

Short Description:

- Easy handling, per rod in the anterior nasal region
- BfArM Test ID: AT282/21
- Sensitivity: 95.37%; Specificity: 99.13%.
- Very easy to use, by swab in the anterior nasal region (Anterior Nasal)
- Test officially approved and suitable for laymen by BfArM
- In single packs

Intended use:

The kit is for the in vitro qualitative determination of SARS-CoV-2 antigen in human anterior nasal specimens. A positive test result indicates that the specimens contain SARS-CoV-2 antigen. A negative test result does not rule out the possibility of infection.

This kit is intended for home use by lay persons outside the laboratory. This includes a person's home or certain non-traditional locations such as offices, sporting events, airports, schools, etc.

Components:

- SARS-CoV-2 Antigen Test Cassette
- Extraction solution in tube
- Disposable Virus Specimen Collection Swabs
- Biohazard sample bags (for safe disposal via household waste)

Features:

- Sensitivity = 95.37%
- Specificity = 99.13%
- Easy to use - Fast and reliable test results in only 15 minutes.
- BfArM listed AT282/21, special approval 5640-S-057/21
- Manufacturer: Beijing Hotgen Biotech Co., Ltd.

Further notes:

- The Hoten Self-Test should be stored at 4 to maximum 30°C.
 - Shelf life is 18 months (date of manufacture and expiry date are indicated on the label).
 - After opening the foil pouch, the test should be used within 30 minutes at room temperature (temperature 10 - 30 °C, humidity $\leq$ 70%).
 - Sensitivity of the test decreases after the first 5 days after the onset of symptoms
- Test can only tell if antigen is present, not in what concentration
- Test result depends on quality of specimen collection, processing, and its transport and storage. Any errors may cause inaccurate test results.
 - People under 18 years of age should be assisted by an adult
 - Intended for in vitro diagnostic use only
 - This product has been approved for the detection of SARS-CoV-2 nucleocapsid protein only, not for other viruses or pathogens.
 - Please be sure to follow the manufacturer's information within the instructions for use and on the product packaging

Unterlagen



Anleitung Hotgen Corona Schnelltest BDA

1 file(s) 1 MB

[Download](#)



Gebrauchsanweisung Hotgen Coronavirus Antigentest

1 file(s) 699 kB

[Download](#)

[Read More](#)

Category: [Corona Self Tests](#)

Tag: [Laientest](#), [Lolly-Test](#), [Schnelltest](#), [Selbsttests](#)



Die fünflagigen FFP2-Maske sind einzeln verpackt und damit ideal zur internen Verteilung geeignet. Farbe weiß und schwarz.

Short Description:

- Protective mask with elastic band and modelable nose clip
- Single packs (reclosable)
- Quality certificate CE 2163 EN 149: 2001 + A1: 2009M
- soft shell
- Color: white & black
- Filtering capacity $\geq 94\%$

Intended use:

This FFP2 mask is the model MA-001 from the company Keholl. The five-layer FFP2 mask is CE certified and complies with the provisions of Regulation (EU) 2016/425 of the European Parliament and of the Council on personal protective equipment and the requirements of the harmonized standard EN 149:2001 + A1:2009. The FFP2 masks reliably protect the wearer and the environment against infection with viruses and against other harmful fine particles.

Features:

- High wearing comfort due to elastic ear straps and a bendable nose clip
- NR - not reusable
- Shape: foldable
- Particle filter efficiency: $\geq 94\%$ (test according to EN 149 FFP2 NR)
- Standards: EN 149:2001 + A1:2009 FFP2 NO.
- Material: 42 % fleece, 29 % hot air cotton, 29 % meltblown fleece
- flexible one size fits all: 3D design, V-shape nose contour and ear loops with elastic waistband ensure a snug fit.
- breathable and skin friendly
- designed for single use and not washable
- GTIN (EAN) barcode for retail sale
- individually packed
- with German labeling
- Brand: Keholl

Further notes:

- This product is for single use only and must be disposed of after use.

A mask which is marked "NR" can only be used once, e.g. for one shift.

[Read More](#)

Category: [Masks](#)

Tag: [Atemschutz](#), [FFP2-Maske](#), [Maske](#)



Der medizinische Mundschutz dient dem Infektionsschutz im Alltag. Der Mundschutz reduziert die Übertragung von Viren von dem Träger an die Luft – er dient somit vorrangig dem Fremdschutz vor Tröpfcheninfektionen

Short description:

- high resistance to liquids and aerosols
- low breathing resistance
- high bacterial filtration
- Typ IIR
- EN 14683:2019 + AC2019
- high bacterial filtration $\geq 98\%$
- 50 pieces per box / 10 pieces per polybag

Intended use:

The surgical mouthguard type IIR from shuwen has an increased resistance to liquids and aerosols as well as a bacteria filtering performance of $\geq 98\%$. With the integrated, long nose clip, the surgical mask can be ideally adapted to the shape of the wearer's face. Fastening with the help of elastic ear loops enables quick donning and doffing. Thanks to the high bacteria filtering performance and resistance to aerosols and liquids, the mouthguard offers particularly high infection protection and is resistant to splashes and spray.

Features:

- Type IIR surgical mouthguard
- Bacterial filtration efficiency $\geq 98\%$
- Color: blue
- 3-layer PP non-woven material
- disposable mask (not reusable)
- incl. adjustable nose clip
- With integrated long nose clip
- Manufacturer: Zhejiang Xichen Medical Technology Co., Ltd.

Further instructions:

- Wash hands and the corresponding part of the face with soap and/or clean with suitable disinfectant before and after removing or otherwise touching the mask
- Wear mask over mouth and nose only
- Remove or change the mask if it becomes damp

[Read More](#)

Category: [Masks](#)

Tag: [Atemmasken](#), [Mask](#), [Medical Mask](#), [MNS](#), [Mund-Nasen-Schutz](#), [OP-Maske](#)



Daddy's Choice Purism FFP2 individually wrapped with very elastic straps and modelable noseband for a comfortable fit.

Short Description:

- Protective mask with elastic band and modelable nose clip
- single packs and boxes of 20
- FFP2 NR (DIN EN 149:2001 + A1:2009)
- CE 0598
- Soft shell mask
- color: white
- Filtering capacity $\geq 94\%$

Intended use:

This FFP2 mask is a Purism brand FFP2 mask. The five-layer FFP2 mask is CE certified and complies with the provisions of Regulation (EU) 2016/425 of the European Parliament and of the Council on personal protective equipment and the requirements of the harmonized standard EN 149:2001 + A1:2009. The FFP2 masks reliably protect the wearer and the environment against infection with viruses and against other harmful fine particles.

Features:

- High wearing comfort due to elastic ear straps and a bendable noseband
- Shape: foldable
- Particle filter efficiency: $\geq 94\%$ (test according to EN 149 FFP2 NR)
- Standards: EN 149:2001 + A1:2009 FFP2 NR
- Flexible one size fits all: 3D design, V-shape nose contour and ear loops with elastic waistband ensure a snug fit
- breathable and skin friendly
- designed for single use and not washable
- Effective filtering system
- Brand: Daddys Choice

Further notes:

- This product is for single use only and must be disposed of after use.
- A mask which is marked "NR" can only be used once, e.g. for one shift.

Unterlagen

[Prüfbericht FFP2-Maske Purism](#)

[EU Zertifikat - FFP2-Maske Purism Daddys Choice](#)

[Read More](#)

Category: [Masks](#)

Tag: [Atemmasken](#), [Atenschutz](#), [FFP2-Maske](#), [Masks](#), [Mund-Nasen-Schutz](#)



Nitril-Handschuhe aus synthetischem Kautschuk sind beständiger gegen Chemikalien und Öle als Latexhandschuhe, undurchlässig gegen Viren und hautverträglicher.

Brief Description:

- Disposable protective gloves
- offers high protection against chemicals, solvents and viruses
- Nitrile rubber (NBR), powder-free, non-sterile
- 100 pieces per box
- CE 2777
- Sizes: S, M, L, XL
- color blue
- ambidextrous

Intended use:

These nitrile gloves are universally applicable. No matter if in medicine, food industry, gastronomy or industry.

Nitrile is a synthetic material that is highly superior to other materials in many areas. Nitrile gloves offer a very good fit and tactile feel, they can be used in many areas. The gloves are used to protect against infections and are highly resistant to chemicals, including acids and alkalis.

Biggest advantage: Nitrile is 100% free of latex proteins and therefore very skin friendly. Nitrile is latex free and therefore allergy free. Ideal for all those who work long and frequently with disposable gloves. Optimal fit and a comfortable wearing feeling are assured.

Features:

- Good skin compatibility - suitable for allergy sufferers
- High flexibility - tear and puncture resistant, but elastic
- Resistant to fats and oils and suitable for foodstuffs
- Sensitivity in the fingertips - good sense of touch
- High abrasion resistance & durability
- Comfortable to wear - good ventilation
- EN 455-1:2020, EN 455-2:2015, EN 455-3:2015 EN ISO 374-1:2016, EN ISO 374-2:2019, EN 16523-1:2015, EN ISO 374-4:2019, 374-5:2016, BS EN ISO 21420, ISO 16604:2004

[Read More](#)

Category: [Gloves](#)

Tag: [Handschuhe](#), [Latexfrei](#), [Nitrilhandschuhe](#), [Puderfrei](#), [Virenschutz](#)



Children FFP2 masks in size S. Perfect for protection against infection due to a small and comfortable fit for children.

Short Description:

- Protective mask with elastic band and modelable nose clip
- 20 pieces per package

- Quality certificate CE 1463 EN 149: 2001 + A1: 2009
- size S - ideal for children
- color: white
- filtering capacity $\geq 94\%$

Intended use:

The five-layer FFP2 mask is CE certified and complies with the provisions of Regulation (EU) 2016/425 of the European Parliament and of the Council on personal protective equipment and the requirements of the harmonized standard EN 149:2001 + A1:2009. The FFP2 masks reliably protect the wearer and the environment against infection with viruses and against other harmful fine particles.

Features:

- High wearing comfort due to elastic ear straps and a flexible nose clip.
- NR - not reusable
- Shape: foldable
- Particle filter efficiency: $\geq 94\%$ (test according to EN 149 FFP2 NR).
- Standards: EN 149:2001 + A1:2009 FFP2 NR
- Flexible one size fits all: 3D design, V-shape nose contour and ear loops with elastic waistband ensure a snug fit.
- breathable and skin friendly
- intended for single use and not washable
- with German label
- Brand: YPHD

Additional Notes:

- This product is for single use only and must be disposed of after use.
- A mask which is marked with "NR" can only be used once, e.g. for one shift.

[Read More](#)

Category: [Masks](#)

